POLYMAILLE®

FLOW

- EN** Instructions for use
- FR** Mode d'emploi
- DE** Gebrauchsanweisung
- ES** Instrucciones de uso
- IT** Istruzioni per l'uso
- PT** Instruções de utilização
- PL** Instrukcja użycia
- TR** Kullanım Kılavuzu

PEROUSE
MEDICAL

Please note instruction for use • Veuillez respecter le mode d'emploi
• Bitte beachten Sie die Gebrauchsanweisung • Por favor, lea las instrucciones de uso • Rispettare le istruzioni per l'uso • Consultar as instruções de utilização • Proszę przestrzegać instrukcji obsługi • Lütfen, kullanım kılavuzuna uyun

Not intended for reuse • Ne pas réutiliser • Nicht zur Wiederverwendung • No reutilizable • Monouso non riutilizzabile • Não reutilizar • Nie używać ponownie • Yeniden kullanmayın.

REF

Order number • Référence produit • Bestell-nummer
• Número de pedido
• Codice d'ordine • Número de referencia • Numer zamówienia • Ürün referansı

Expiry date • Date de péremption • Verwendbar bis • Fecha de caducidad • Data di scadenza • Prazo de validade • Data ważności • Kullanım süresi dolumu

Date of manufacture • Date de fabrication • Herstellungsdatum • Fecha de fabricación • Data di produzione • Data de fabrico • Data produkcji • Üretim tarihi

STERILE|EO

Sterilized using ethylene oxide
• Stérilisé à l'oxyde d'éthylène • Sterilisation mit Ethylenoxid • Esterilización con óxido de etileno • Sterilizzazione con ossido di etilene • Esterilizado com óxido de etileno • Steryliizowane tlenkiem etylenu
• Etilen oksitile sterilize edilmişdir

English **P. 4**

Français **P. 12**

Deutsch..... **P. 21**

Español **P. 30**

Italiano..... **P. 38**

Português **P. 47**

Polski **P. 56**

Türkçe **P. 65**

POLYMAILLE® FLOW

ePTFE Vascular Prosthesis

Product description

POLYMAILLE® FLOW is a vascular prosthetic graft made of expanded polytetrafluoroethylene (ePTFE). A thin ePTFE wrap is wound around the single-layer body of the graft.

The graft is available in two wall thicknesses, „thin wall“ and „standard wall“. It is also available with a removable external helical Polytetrafluoroethylene (PTFE) reinforcement for sections in which the risk of compression and kinking is high (reconstructions that traverse a joint or extraanatomical reconstructions).

Four stickers for patient documentation are enclosed in the packaging.

Indications

The POLYMAILLE® FLOW is indicated in cases of arterial vasoreconstruction, primarily in the peripheral vascular region. It can also be used in extraanatomical reconstructions - femorofemoral and axillofemoral.

Standardwall POLYMAILLE® FLOW grafts are also indicated for use as arteriovenous shunt prostheses in haemodialysis.

Contraindications

The available clinical data are insufficient to support definitive conclusions concerning the use of ePTFE vascular grafts in the following areas of application:

Reconstruction

- of the pulmonary arteries
- of the aorta ascendens
- of the coronary arteries
- of the carotid arteries (A. carotis communis, externa and interna)
- of the cerebral arteries
- of the arteria brachialis and cephalica
- of the cardiac and pulmonary veins and
- of the vena cava inferior and/or superior

The POLYMAILLE® FLOW should also not be implanted in patients who show an immunological hyperreaction to ePTFE.

Precautions and warnings

Failure to comply with these precautions and warnings may result in serious medical complications.

- The POLYMAILLE® FLOW ePTFE vascular prosthetic graft should only be implanted by persons with special training in the surgical implantation of the material ePTFE.
- Do not use this product after the expiry date. Please ask your supplier in such cases.
- The graft must not be exposed to strong radiation. ePTFE decomposes when exposed to high doses of radiation.
- Do not allow the graft to come into contact with reagents containing oil, alcohol or aqueous solutions. Do not attempt to preclot the graft.

- Keep processing and manipulation of the graft to a minimum, since otherwise blood, serum or plasma could be pressed through its wall.
- Useatraumatic clamps with soft, cushioned ends.
- Always use non-absorbable suture material. The suture seam should not measure less than 2 mm on vessels.
- Thin-walled ePTFE vascular prostheses are generally not recommended for extraanatomical, axillofemoral vessel prosthesis, since the anastomosis may rupture under heavy stress loads.
- POLYMAILLE® FLOW grafts with helical reinforcement over their entire length are not indicated for use as arteriovenous shunt prostheses. In the graft models with central helical reinforcement, blood access is not permissible in the area of the reinforcement. There is an increased risk of embolism and/or pseudoaneurysms in both cases.
- Do not use this prosthetic graft or parts of it as a patch. The product design is not suitable or beneficial for applications other than those indicated and undesirable events may result from such misuse.
- POLYMAILLE® FLOW is a single-use product that has been sterilized with EtO. It may not be resterilized or reused. Risks arising from resterilization and reuse include infection, loss of function, hemolysis, allergic reaction in the patient. The product must not be used, if the package is damaged or if the package has been opened outside of the sterile area.

Adverse reactions

As with any surgical implant, prosthetic infection necessitating explantation is a possibility. Maintenance of stringent aseptic surgical conditions is therefore of the utmost importance.

Other potential complications include thrombosis, suture haemorrhaging along the anastomosis, suture rupture, elongation, seroma formation, development of pseudoaneurysms and haematomas.

Risks involved in use as an arteriovenous shunt graft: steal syndrome, skin erosion.

Implantation procedure

Artificial vascular prostheses must always be implanted by surgeons trained in and familiar with the specific techniques of such operations. Please also observe the following points when implanting the POLYMAILLE® FLOW ePTFE vascular prosthetic graft:

Preparation:

- When opening the packaging, please be aware that it is a double-blister. The inside blister only may be brought into the sterile field.
- Touch the graft with sterile, unpowdered gloves only. Use no sharp or pointed instruments and atraumatic vascular clamps only.
- When selecting the diameter, make sure the dimension matches the lumen of the original physiological blood vessel as closely as possible.

- Before cutting the prosthesis to size, first remove a suitable length of the spiral reinforcement along the outside from each end as follows: Take hold of the end of the spiral with the fingers or an atraumatic instrument, then slowly and carefully unwind the spiral reinforcement, holding it at a 90° angle to the prosthesis.
- Cut the graft to size use a sharp scissor or scalpel. Since ePTFE is not elastic in the longitudinal direction, make sure the cut length is correct, taking into account the patient's position, body weight and ensuring complete freedom of movement of the extremity or extremities. If the graft is cut off too short, the suture holes may be stretched, leading to suture haemorrhaging. Cutting it too long may result in kinks.

Placement using tunneler:

- Helix-reinforced grafts require a tunneler with a bullet tip 1 - 2 mm larger than the inside diameter of the graft to facilitate unobstructed passage and placement. Before pulling the graft through the tunnel, replace the bullet tip with one that corresponds exactly to the inside diameter of the graft.
- When using grafts without external reinforcement, the diameter of the tunneler should be the same as the inside diameter of the graft to be implanted.

Anastomosis:

- Use only atraumatic clamps with soft, cushioned ends.
- For suturing use a non-cutting round needle of suitable size and a corresponding non-absorbable monofilament suture material.

- The suture seam along the anastomotic graft tissue should have the same width on both sides. This seam margin should not be less than 2 mm in reconstructions of vessels.
- Avoid excessive tension on the anastomosis as well as inappropriate suture gaps.
- For axillofemoral bypasses, the recommended location for the proximal anastomosis is as far medial as possible in the first third of the a. axillaris. Perform the anastomosis procedure with the patient's arm outstretched.
- Minimize the risk of suture haemorrhaging by following the curve of the needle, holding the suture thread at a 90° angle as far as possible and never pulling hard on the suture seam.

Thrombectomy:

- POLYMAILLE® FLOW ePTFE vascular prosthetic grafts are not radially elastic. For postoperative placement of balloon catheters and/or angioplasty catheters in the graft lumen you must therefore make sure the outside diameter of the inflated balloon is the same as the inside diameter of the graft. Excessive inflation and/or use of the wrong size may result in dilatation of and/or damage to the prosthesis.

POLYMAILLE® FLOW ePTFE vascular prosthetic graft used as an arteriovenous shunt graft in haemodialysis:

- To reduce haematoma formation, blood access (catheterization) through the shunt should not be carried out until about 2 weeks after the implantation procedure.

- Make sure the dialysis cannula is inserted at a 45° angle to the graft, then position it parallel to the prosthesis.
- Do not cannulate the same insertion point twice, since repeated puncturing in the same area may result in formation of haematomas and pseudoaneurysms.

Storage and packaging

The POLYMAILLE® FLOW ePTFE vascular prosthetic graft is an EO-sterilized, single-use product that is supplied in a sterile package. It will remain sterile as long as the package is neither opened nor damaged.

Store in a cool, clean, dry place.

Replacement warranty

Please ship the product with all packaging and accessories back to your supplier without delay in case of accidental opening of or damage to the packaging or if the packaging, or indeed the product, has been damaged in transit. The supplier will see to it that the items are replaced immediately.

WARRANTY AND LIABILITY LIMITATIONS

PEROUSE MEDICAL warrants that this product conforms to the legal requirements, in particular those set forth in the EU Guidelines for Medical Products, 93/42/EWG. Therefore, in case this product exhibits faults that can be directly attributed to material or processing errors that occurred during its manufacture or packaging („faults“), and said faults appeared prior to the expiration date (UBD), PEROUSE MEDICAL will, at its discretion, replace the product

- a) with one of equal construction, or;
- b) with one exhibiting the identical function.

All the following conditions must be met in order to make claims within the context of this warranty:

- a) the product must not have been employed after its UBD.

- b) the fault did not initially occur after the UBD.
- c) if the product in question is a disposable product or is a product which PEROUSE MEDICAL designates as being for single use, the fault must be detected at the first use.
- d) the fault was immediately reported to PEROUSE MEDICAL or its representative.
- e) the faulty product was forwarded to PEROUSE MEDICAL or its representative immediately after the fault was reported.

This warranty in no way limits and applies in conjunction with any and all regulations governing warranties on the part of the seller and the product liability laws.

For the remainder, this warranty applies only to the extent defined above. Additional claims arising from this warranty are excluded. This warranty as well as the cited legal claims apply in place of, and to the exclusion of, all other warranties, either verbal or written, or express or implicit on the part of PEROUSE MEDICAL. None of our representatives or employees are empowered to reach differing agreements.

Due to numerous factors such as shipping, storage, and handling on the part of the user, all of which lie outside the control of PEROUSE MEDICAL, PEROUSE MEDICAL is not liable, in particular, for damages of any type resulting directly or indirectly from the product or its employment unless said damages can be shown to be the result of a failure on the part of the product.

In particular, PEROUSE MEDICAL assumes no liability for damages resulting from resterilization or reuse of products approved for single use only, or which PEROUSE MEDICAL designates as being for single use.

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