# POLYSITE®

- EN** Implantable ports
- FR** Chambres implantables
- DE** Implantierbare Portsysteme
- ES** Cámaras implantables
- PT** Câmaras implantáveis
- NL** Implanterbare poorten
- TR** İmplemente edilebilir portlar
- AR** حجيرات قابلة للغرس

*INSTRUCTION AND INDICATIONS*

*INSTRUCTIONS ET INDICATIONS*

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## INSTRUCTIONS AND INDICATIONS

### I/ DESCRIPTION / INDICATIONS

POLYSITE® implantable access ports allow repeated and prolonged access to the vascular system (intravenous or intra-arterial), peritoneal cavity or intra-vertebral space (intrathecal and peridural space).

- POLYSITE® venous access ports are used to administer chemotherapy, antibiotics and antiviral drugs. They can also be used for parenteral nutrition, collection of blood samples and transfusion of blood or blood products. Some POLYSITE® references can be used for high flow rate injection of contrast media during CT procedures (compatible with high pressure injectors, cf Chapter VII/ B/). The POLYSITE® ECHO placement kits are used in the same indications and consist of instrumentation suitable for placement of implantable ports by echo-guided puncture (see chapter III).
- POLYSITE® arterial and peritoneal access ports are used for local or regional chemotherapy. Blood samples cannot be taken via arterial access ports, as they are equipped with an anti-reflux valve.
- POLYSITE® spinal access ports are used to administer analgesics and cytotoxic drugs.

POLYSITE® implantable access ports are composed of a radiopaque port (or housing) connected to a catheter and maintained by a connection ring. They are supplied either non-assembled or pre-assembled. The port is accessed percutaneously by using a Huber needle.

The materials of which the POLYSITE® implantable ports are made can be used with almost any injectable solution (c.f. Chapter V/ A/) and are compatible with CT and MRI procedures.

### II/ CONTRAINDICATIONS

Inserting of a POLYSITE® implantable access port is contraindicated in the following situations:

- Known or suspected history of immunoallergic or hypersensitivity reactions to materials contained in the device.
- Presence of local or systemic infection, bacteraemia or septicaemia.
- When the substances used for treatment are incompatible with one of the materials of the device.

### III/ CONTENT

- Each POLYSITE® pack contains:
  - An acetal resin [Polyoxymethylene or POM] or VOLTALEF® [PCTFE] or titanium and POM port with a silicone septum.
  - One silicone, or XRO polyurethane catheter.
  - Two connection rings (only one for pre-assembled models).
  - A straight Huber needle.
  - A vein holder for POLYSITE® venous access ports.
  - A Tuohy needle and a PTFE-coated metal guide-wire (mounted inside the catheter) for POLYSITE® spinal access ports.
  - One catheter flush for POLYSITE® venous access ports, not pre-connected.
- Each POLYSITE® kit also contains:
  - An introducer (dilator with peelable or non-peelable sheath).
  - A J-shaped guidewire.
  - A puncture needle.
  - A tunneler.
  - One syringe.

- Each POLYSITE® ECHO placement kit contains:
  - An acetal resin [POM] or VOLTALEF® [PCTFE] or titanium and POM port with a silicone septum.
  - One silicone, or XRO polyurethane catheter.
  - Two connection rings (only one for pre-assembled models).
  - A straight Huber needle.
  - A vein holder for POLYSITE® venous access ports.
  - One catheter flush for POLYSITE® venous access ports, not pre-connected.
  - An introducer (dilator with peelable or non-peelable sheath).
  - A J-shaped guidewire.
  - An echogenic puncture needle.
  - A tunneler.
  - A 10 ml syringe.
  - A sterile drape.
  - A Raulerson device introducer.
  - A non-latex echo-guidance probe sheath.
  - A sterile ultrasound gel sachet.
  - Two supporting elastic straps for the protective cover.

#### **IV/ POSSIBLE COMPLICATIONS**

The complications most frequently described after insertion of an implantable access port include:

- Gas embolism.
- Pneumothorax, haemothorax, hydrothorax.
- Haemorrhage, haematoma.
- Brachial plexus injury.
- Cardiac arrhythmia, tamponade.
- Thoracic duct injury.
- Arteriovenous fistula.
- Occlusion and/or rupture of the catheter due to compression between the clavicle and the first rib.
- Disconnection, rupture, embolisation of the catheter.
- Occlusion of the catheter or port.
- Thromboembolism/phlebitis, thrombosis.
- Drug extravasation.
- Cutaneous erosion.
- Local or systemic infection: bacteraemia and/or septicaemia.
- Formation of a fibrin cuff.
- Migration of the implant, inversion of the port.
- Inflammation, skin necrosis over the implantation site.
- Abdominal hernia, peritonitis, peritoneal leak (peritoneal sites).
- Intolerance reactions to the device.

#### **V/ WARNINGS AND PRECAUTIONS**

The warnings and precautions described below must be rigorously observed to avoid any incidents likely to cause device dysfunction and a risk to the patient.

##### **A/ GENERAL WARNINGS**

- ➔ Disposable medical device. **DO NOT REUSE.** The POLYSITE® implantable ports must never be re-implanted. As they are contaminated by blood a used implantable port represents a biological risk (risk of infection) In addition, PEROUSE MEDICAL cannot guarantee the correcting functioning of the device after it has been removed (system patency, septum watertightness, catheter integrity, etc.)
- ➔ Sterilized with ethylene oxide. **DO NOT RESTERILIZE.**

- **Do not use the device if the packaging is damaged, opened or after the expiry date.**
- Store protected from light, cold, heat ( $10^{\circ}\text{C} < T < 35^{\circ}\text{C}$ ) and humidity (RH < 80%). Under these conditions, the product can be used until the expiry date indicated on the packaging.
- Observe a strict aseptic technique during implantation and use of the device.
- Do not use medicinal products or solvents which are incompatible with the catheter material. Note: do not use alcohol with polyurethane catheters.

## **B/ INSERTION**

- Only qualified doctors can insert an implantable access port.
- Before insertion, purge the device of any air and fill with physiological saline or heparinised solution.
- When using the POLYSITE® ECHO placement kit, **PURGE THE RAULERSON DEVICE INTRODUCER WITH PHYSIOLOGICAL SALINE OR HEPARINISED SOLUTION BEFORE PUNCTURING THE VEIN to prevent air from entering the vein (risk of gas embolism).**
- Take care during insertion of the device to avoid any mechanical alteration of the catheter. **Only clamp the catheter withatraumatic clamps or forceps. Do not use the catheter if it is damaged or if it leaks.**
- **The catheter and port must always be assembled manually.** Strictly follow the connection technique described under VI/ INSERTION TECHNIQUES to avoid damaging the catheter and to ensure correct connection.

### ***Subclavian catheter insertion***

**A subclavian catheter can be compressed between the clavicle and the first rib (pinch-off phenomenon) with a risk of rupture and embolisation of the catheter.**

Excessive and/or repeated movements of the upper limb increase this risk.

- A number of signs are suggestive of pinching of the catheter, including the need to mobilize or raise the patient's arm to perform stable infusion.
- To avoid this phenomenon, catheterize the subclavian vein away from the costoclavicular angle.
- Take an x-ray to make sure that the catheter is not pinched.
- Remove the catheter in the presence of signs of pinching.

## **C/ USE AND MAINTENANCE**

- **Lack of maintenance can lead to occlusion of the device.**
- **Only use Huber needles or POLYPERF® or POLYPERF Safe® connection lines.** Conventional hypodermic needles should not be used, as they can alter the integrity of the septum.
- **Do not reuse** Huber needles or POLYPERF® connection lines. Deformation of the bevel of the needle can tear or damage the septum of the port.
- Only use large gauge Huber needles (19 and 20 G) for blood transfusions, parenteral nutrition and injection of contrast media.
- Select the length of the needle by taking into account the depth of the port, the thickness of the tissues and the thickness of the dressings underneath the curvature of the needle. Needles that are too short or too long may be displaced by the patient's movements with a risk of drug extravasation.
- **NEVER TRY TO UNBLOCK THE DEVICE BY PRESSURE AND NEVER USE A SYRINGE SMALLER THAN 10 ml to avoid applying excessive pressure on the device.** Excessive pressure can damage the catheter or lead to disconnection of the port.
- When the implantable port is used to inject contrast media, first confirm that it is completely patent (by obtaining blood reflux) and follow the instructions listed in chapter VII/ B/.
- Always ensure that the POLYSITE® implantable port reference has been validated for use in injecting contrast media (cf table 1) before using with a high pressure injector. If the POLYSITE® implantable port reference has not been validated for this use, high pressure

injection of contrast media risks damaging the system (catheter rupture, disconnection of port from catheter and cause complications to the patient (extravasation of contrast media, migration of the catheter).

- Do not use Huber needles with incompatible tubing which has not been validated for use with a high pressure system.
- Do not administer the injection with a high pressure injector if the port is at risk of compression or costo-clavicular compression.
- Do not administer high pressure injections through catheters over 25 cm long because of a risk of significantly reducing the flow rate.

#### D/ DURATION OF IMPLANTATION AND REMOVAL OF THE IMPLANTABLE PORT

- The implantable catheter port must be removed at the end of treatment.
- The lifespan of the system depends on the duration of treatment, tailored to each disease and the functional state of the device, particularly the catheter. The most common complication of ageing is fracture.
- **Do not leave the device in place when it is no longer necessary.**
- **Remove the catheter before removing the port.** Do not pull on the port while the catheter is still inserted, as this can lead to disconnection of the device and embolisation of the catheter.

### VI/ INSERTION TECHNIQUES

#### A/ IMPLANTATION OF POLYSITE® VENOUS ACCESS PORTS

The doctor performing insertion must determine the most appropriate method of introduction for the patient concerned, after having evaluated the risks and advantages of each method.

The risk of gas embolism can be avoided by placing the patient in the Trendelenburg position (head-down position).

##### Surgical method

1. Purge the port and the catheter with physiological saline or heparinised solution.
2. Perform venous cutdown.
3. Introduce the catheter into the vein using a vein holder. Suture of the catheter to the vein is useless or even dangerous (risk of venous stripping during removal).
4. Advance the catheter as far as the junction of the superior vena cava and right atrium.
5. Check the position of the distal extremity of the catheter on an x-ray.
6. Prepare a subcutaneous site so that the port is situated about 1 cm under the surface of the skin.
7. Perform trial insertion of the port to check that the dimensions of the site are suitable and that the port is not situated underneath the incision.
8. Tunnel the catheter as far as the port site.
9. Cut the proximal extremity of the catheter to the required length, making sure that there is sufficient slack to allow body movements and connection to the port.
10. Thread the connection ring over the catheter.
11. Mount the catheter onto the port to cover the conical tip (see diagrams below).
12. Push the ring (which advances the catheter) until it locks. The tip of the catheter must not be compressed over the port.

POLYSITE® 2000, 3000 & 4000

POLYSITE® 1000

13. Check good connection of the catheter by applying slight traction.
14. Insert the port into the port site and check for the absence of kinking of the catheter.
15. Check the patency of the device.
16. Anchor the port in its site by using permanent sutures, and then close the incisions.
17. Rinse and heparinise the device.

#### Percutaneous method

1. Purge the port and the catheter with physiological saline or heparinised solution.
2. Puncture the vein with the puncture needle connected to a syringe.
3. As soon as blood reflux is observed, remove the syringe and leave the needle in the vein.
4. Insert the metal guide-wire into the vein. Remove the needle.
5. Enlarge the cutaneous puncture site to allow passage of the introducer.
6. Thread the introducer over the guide and advance it by exerting a rotation movement.
7. Remove the metal guide-wire and the dilator of the introducer. Leave the sheath in place.
8. Insert the catheter into the sheath.
9. Advance the catheter as far as the junction of the superior vena cava and right atrium.
10. Check the position of the distal extremity of the catheter on an x-ray.
11. Remove the sheath.
12. Follow points 6 to 17 of the surgical method described above.

#### Percutaneous method using the echo-guided puncture technique

1. Purge the port and the catheter with physiological saline or heparinised solution.
2. Apply sterile gel to the end of the probe, then fit the protective sheath over it. The sheath is held on the ultrasound probe by two rubber bands supplied in the kit.
3. Apply a layer of sterile gel over the protective sheath at the level of the probe end.
4. Apply the ultrasound probe to the skin surface over the chosen venous access point until this can be viewed on the ultrasound monitor.
5. Connect the Raulerson device to the puncture needle.
6. **BEFORE PUNCTURING THE VEIN, to prevent air from entering the vein (risk of gas embolism), INSERT THE PUNCTURE NEEDLE CONNECTED TO THE RAULERSON DEVICE INTO A BOTTLE PHYSIOLOGICAL SALINE (upside down), AND SQUEEZE SEVERAL TIMES THE RESERVOIR OF THE RAULERSON DEVICE UNTIL COMPLETELY FREE OF AIR.**
7. Puncture the vein with the puncture needle fitted to the Raulerson device kept squeezed.
8. As soon as the point of the puncture needle is in the venous access, the reservoir of the Raulerson device fills with blood (venous reflux). Proper positioning of the needle point in the venous access can be checked on the ultrasound monitor.
9. Do not disconnect the Raulerson device from the puncture needle.
10. Insert the guidewire into the vein through the Raulerson guidewire introducer and puncture needle. When inserting the guidewire through the lumen of the Raulerson device, the user can feel some resistance corresponding to the passage of the guidewire through the check valve.
11. Remove the puncture needle and Raulerson device.
12. Follow points 1 to 12 of the percutaneous method described above.

#### Specificities concerning insertion of pre-assembled access devices

1. Before inserting the catheter into the vein:
  - Prepare the subcutaneous port site.
  - Tunnel the catheter from the port site.
  - Cut the catheter to the required length at its distal extremity.
2. A peelable introducer must always be used for percutaneous insertion.

## **B/ INSERTION OF POLYSITE® ARTERIAL ACCESS PORTS**


1. Perform arteriography before insertion.
2. Purge the port and the catheter with heparinised solution.
3. Perform arterial cutdown.
4. Introduce the catheter into the artery via its distal extremity.
5. Check the correct position of the catheter in the artery and permanently ligate using fixation rings.
6. Check organ perfusion.
7. Follow points 6 to 17 of the surgical method described for POLYSITE® venous access ports.

## **C/ INSERTION OF POLYSITE® PERITONEAL ACCESS PORTS**

1. Purge the port and the catheter with physiological saline.
2. Perform an incision through the abdominal wall as far as the peritoneum.
3. Introduce the distal extremity of the catheter into the abdominal cavity as far as the desired infusion site.
4. Check the position of the catheter.  
**To avoid infusion into the subcutaneous space, make sure that all side orifices of the catheter are situated in the peritoneal cavity.**
5. Follow points 6 to 16 of the surgical method described for POLYSITE® venous access ports.
6. Rinse the device.

## **D/ INSERTION OF POLYSITE® SPINAL ACCESS PORTS**

1. Purge the port and the catheter with physiological saline.
2. Perform an incision between the spinous process of L2 and L3 or L3 and L4.
3. Introduce the Tuohy needle into the centre of the incision.
4. Introduce the catheter, over its metal guide-wire, through the Tuohy needle.
5. Advance the catheter with image intensifier guidance to the desired position.
6. Remove the metal guide-wire from the catheter.
7. Remove the Tuohy needle while maintaining the catheter in place.
8. Prepare a subcutaneous site so that the port is situated about 1 cm under the surface of the skin.
9. Perform trial insertion of the port to check that the dimensions of the site are suitable and that the port is not situated underneath the incision.
10. Tunnel the catheter as far as the port site.
11. Cut the proximal extremity of the catheter to the required length with a long bevelled cut, making sure that there is sufficient slack to allow body movements and connection to the port.
12. Thread the connection ring over the catheter. The catheter must sufficiently overlap the ring for easy assembly on the connector.
13. Mount the catheter onto the neck of the connector (see diagram below).
14. Push the ring (which advances the catheter) until it locks. The tip of the catheter must not be compressed over the port.

15. Insert the port into the port site and check for the absence of kinking of the catheter.
16. Check the patency of the device.
17. Anchor the port in its site by using permanent sutures, and then close the incisions.

## VII/ USE AND MAINTENANCE

### A/ GENERAL INFORMATION

1. Before use, inspect and perform asepsis of the injection site.
2. Locate the port by palpation.
3. Insert the needle perpendicularly to the skin and advance as far as the base of the port.
4. Check the patency of the device and correct positioning of the needle by obtaining blood reflux (POLYSITE® venous access ports only) and by injection of physiological saline without excessive pressure or local extravasation.

**When in doubt about the integrity of the device, defer chemotherapy and perform opacification of the device with contrast agent.**

5. Inject or infuse the drug solution or take the blood sample (eliminating the first few ml of blood taken).

If several products must be administered, rinse with physiological saline after each injection.

**At the slightest sign of extravasation, stop injection immediately and take appropriate measures.**

6. Rinse the device with physiological saline after each use (this operation is unnecessary for spinal and peritoneal access ports).
7. After taking blood samples rinse with physiological saline and if necessary flush with heparin (see the protocol of the establishment). No maintenance is required for peritoneal or spinal access ports.
8. Remove the needle by maintaining positive pressure to avoid reflux of blood inside the catheter.

### B/ INJECTION OF CONTRAST MEDIA DURING CT SCAN PROCEDURES

1. Follow stages 1 to 4 described above (USE AND MAINTENANCE, A/ GENERAL INFORMATION).
2. Always ensure that the POLYSITE® implantable port reference is present in table 1 before using a high pressure injector.
3. Use a Huber needle without tubing or a Huber needle with tubing validated for high pressure injection.
4. Select the appropriate Huber needle gauge to the viscosity of the contrast medium used, the and the implantable port reference following the recommendations in table 1.
5. Set the maximum injector pressure to 22.4bars or 325 psi in order to guarantee the reliability of the system.
6. Only use 19G or 20G needles for POLYSITE® Standard (series 4000) and Mini (series 3000) and 20G or 22G for POLYSITE® Micro (series 2000), in accordance with table 1.
7. Do not inject more than 3 ml/s when using a 22G Huber needle.
8. Always ensure that the injection circuit is fully operational (Huber needle and implantable port) by obtaining reflux of blood and injecting 10-20 ml of physiological saline without difficulty.
9. **Warm the contrast medium to 37°C (100°F) before use.** Failure to observe this measure would reduce the flow rate by approximately 50% and cause malfunction of the port or injection system.
10. **Always ensure that the catheter is equal or less than 25cm long** (do not administer this type of injection with a femoral catheter). Flow rate may be significantly reduced if a longer catheter is used.
11. **Never exceed the maximum recommended flow rate for a given port.**
12. Never inject contrast medium with a viscosity greater than that shown in table 1.
13. Rinse the implantable port with 10 to 20 ml of 0.9% NaCl before and after use followed by usual rinsing procedures.

**Failure to follow these recommendations risks causing failure of the system through excess pressure or obstruction (c.f. Chapter V/ WARNINGS AND PRECAUTIONS).**

Table 1 shows the maximum flow rates which we recommend for injecting contrast medium and the precise conditions to follow for this use.

## VIII/ TRACEABILITY

Self-adhesive labels and a patient card are provided in each packaging to allow the required follow-up of implanted POLYSITE® access ports. One label must be pasted onto the patient's medical record and another onto the patient card. A medical device vigilance report form and a surveillance notebook to be given to the patient are also enclosed.

**TABLE 1**

Maximum flow rates which we recommend for injecting contrast medium, the following conditions must be observed:

	POLYSITE® reference	Maximum recommended flow rate (mL/s) with 25 cm catheter	Maximum recommended pressure (CT scan function)	Viscosity (cP) / maximum recommended iodine concentration (mgI/mL) in contrast media	Recommended diameter of angled needle (Gauge)			
Micro Series	2005	1	325 psi or 22,4 bars	6cP i.e. 300 mgI/mL (e.g.: Xenetix® 300, pre-warmed to 37°C)*	20G and 22G			
	2005 ISP, 2005 ECHO							
	2105, 2105 ISP, 2105 ECHO	3						
	2016, 2016 ISP, 2016 ECHO							
	2116	3						
	2116 ISP, 2116 ECHO							
	2007							
	2007 ISP, 2007 ECHO	5						
Mini Series	2107							
	2107 ISP, 2107 ECHO							
	3007	3	10cP i.e. 350 mgI/mL (e.g.: Xenetix® 350, pre-warmed to 37°C)*	19G and 20G				
	3007 ISP, 3007 IS, 3007 ECHO							
	3107							
	3107 ISP, 3107 IS, 3107 ECHO							
	3017							
	3017 ISP, 3017 ECHO	5						
	3117							
	3117 ISP, 3117 ECHO							
Standard Series	3008	5						
	3008 ISP, 3008 IS, 3008 ECHO							
	3108							
	3108 ISP, 3108 IS, 3108 ECHO							
	4017	5	10cP i.e. 350 mgI/mL (e.g.: Xenetix® 350, pre-warmed to 37°C)*	19G and 20G				
	4017 ISP, 4017 ECHO							
	4117							
	4117 ISP, 4117 ECHO							
	4008	5						
	4008 ISP, 4008 ECHO							
	4108							
	4108 ISP, 4108 ECHO							
	4019	5						
	4019 ISP, 4019 ECHO							
	4119							
	4119 ISP, 4119 ECHO							
	40010	5						
	40010 ISP, 40010 ECHO							
	41010							
	41010 ISP, 41010 ECHO							

\* Xenetix® is a registered trade mark of Laboratoires Guerbet.

**TABLE 2**

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POLYSITE® is a registered trademark of PEROUSE MEDICAL