

Huber Needles and Polyperf®

COMPLETE RANGE GAUGES SIZES

For the puncture of venous, intra-rachidian, arterial and intra-peritoneal implantable ports

Huber needles

References	Description			Usable length (mm)	Recommended field of use	Packaging
522506	Curved	6/10°	23G	20	Pediatric chemotherapy	Box of 50
522007		7/10°	22G	20	Chemotherapy, Blood reflux, Hyper-hydration	
522507		7/10°	22G	25		
523807		7/10°	22G	35		
522509		9/10°	20G	20	Hyper-hydration, Chemotherapy, Packed red blood cells, Blood collection, Nutritional fluids delivery	
523809	20G		35			
522511	11/10°	19G	20	Heavy flows, Packed red blood cells, Blood collection, Transfusion, Nutritional fluids delivery		
523811		19G	35			
512506	Straight	6/10°	23G	25	Pediatric chemotherapy	Box of 50
512507		7/10°	22G	30	Chemotherapy, Blood reflux, Hyper-hydration	
513807		7/10°	22G	38		
512509		9/10°	20G	25	Hyper-hydration, Chemotherapy, Packed red blood cells, Blood collection, Nutritional fluids delivery	
512511		11/10°	19G	25		
513811	19G		38			

Ethylene oxide sterilization

Large winglets : - good grip
- good stability

TUBULAR - Non DEHP (without phthalates)

No interaction with anti-cancerous agents

Polyperf® without lateral injection site - Huber needles

References	Description			Usable length (mm)	Tubular (cm)	Packaging
581707	Curved needle with luer lock and clamp	7/10° mm	22G	17	25	Box of 10
582007		7/10° mm	22G	20	25	
582507		7/10° mm	22G	25	25	
583007		7/10° mm	22G	30	25	
583507		7/10° mm	22G	35	25	
581709		9/10° mm	20G	17	25	
582009		9/10° mm	20G	20	25	
582509		9/10° mm	20G	25	25	
583009		9/10° mm	20G	30	25	
583509		9/10° mm	20G	35	25	
581711		11/10° mm	19G	17	25	
582011		11/10° mm	19G	20	25	
582511		11/10° mm	19G	25	25	
583011		11/10° mm	19G	30	25	
583511		11/10° mm	19G	35	25	


Ethylene oxide sterilization

Polyperf® with lateral injection site - Huber needles

References	Description			Usable length (mm)	Tubular (cm)	Packaging
592007	Curved needle with luer lock and clamp with lateral injection site	7/10° mm	22G	20	25	Box of 10
592507		7/10° mm	22G	25	25	
592009		9/10° mm	20G	20	25	
592509		9/10° mm	20G	25	25	
592011		11/10° mm	19G	20	25	
592511		11/10° mm	19G	25	25	

Ethylene oxide sterilization

CE 0459



Winglets colour coded :
- easy identification
of the gauge

