


Cardiovascular

▶ POLYPATCH®

Thin wall knitted polyester
collagen coated

POLYPATCH

- ▶ Easy handling
- ▶ Anatomical shape
- ▶ Resist suture line bleeding



FORMALDEHYDE FREE

PEROUSE
MEDICAL

Up to excellence

POLYPATCH

INDICATIONS

- Carotid angioplasty
- Femoral angioplasty

PROPERTIES

- Permeability : <math><10 \text{ ml/cm}^2/\text{min}</math> under 120 mm Hg
- Thickness : 0,45 mm
- Suture transversal resistance : $\geq 17 \text{ N}$
- Coated bovin collagen
- Sterilized by β irradiation

REFERENCES

References	Designation	Dimensions
PA06075	Anatomical	6 x 75 mm
PA08075	Anatomical	8 x 75 mm
PA10075	Anatomical	10 x 75 mm
PR06075	Rectangular	6 x 75 mm
PR08075	Rectangular	8 x 75 mm
PR10075	Rectangular	10 x 75 mm
PR14075	Rectangular	14 x 75 mm
PR20120	Rectangular	20 x 120 mm

Individually packaged

This device is CE marked.