


Interventional
imaging

▶ DOLPHIN

The inflation device matching
every catheter balloon

- 30 atm pressure gauge
- 30 cc syringe
- Ergonomic handle



PEROUSE
MEDICAL

Up to excellence

▶ DOLPHIN



● PRESSURE GAUGE

- 30 atm
- ▶ reliable in any procedures
- fluorescent indicator
- ▶ good pressure reading

● BODY AND BUTTON


- 30 cc syringe
- ▶ small to big balloons (PTA & PTCA)
- square shaped
- ▶ stable
- automatic cam locking system
- ▶ easy handling
- clear material
- ▶ easy air bubble detection and removal

● HANDLE

- T-shape handle
- ▶ high torque
- ▶ easy deflation

References	Description	Accessories	Packaging
0185ND	30 atm inflation device		6-unit box
0185NA	30 atm inflation device	trefoil 3-way stopcock	6-unit box
0185QL	30 atm inflation device + PTCA set 3 items	trefoil 3-way stopcock myshell lite Y connector 9F 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer alligatork torque device	6-unit box

Please contact **PEROUSE MEDICAL** for non sterile and OEM products



PEROUSE
MEDICAL

135 Route Neuve – 69540 Irigny – France
Tél: 33 (0) 4 72 39 74 10 - Fax: 33 (0) 4 78 51 89 67

www.perousemedical.com