

pelicangio



## A RANGE OF SYRINGES for contrast media injectors

Each case is different ...

- A wide size range including 100, 130, 150 and 200 ml allows **pelicangio** syringes to meet all your needs.
- **pelicangio** syringes match the main injectors available on the market: Medrad, Covidien, Nemoto and Medtron for vascular, CT and MRT.

Operating safety is your priority ...

- **pelicangio** syringes provide a protection against air bubble injection with a dual lip seal secured to the plunger and preventing air aspiration.
- The high quality and very clear material of **pelicangio** syringes makes air bubbles detection easier.
- A protective cap is fitted to prevent contamination until the syringe is being used.



ERGONOMIC APPROACH  
ENHANCED SAFETY  
EASY TO USE

**PEROUSE**  
MEDICAL



## TECHNICAL FEATURES

- Fast and easy connection for every standard Luer lock
- High injection pressure resistance (up to 1200 PSI or 83 bars)

REFERENCES	VOLUME	APPLICATION	INJECTORS	ACCESSORIES	PACKAGING
<b>Medrad injectors</b>					
0137ND	200 ml	CT	CT 200/300; MCT/MCT + ; OP-100; VCT610; ECT710	Fill Tube	50-unit box
0137SD	60 ml	Angio.	Mark V, V+ & Provis	Fill Tube	50-unit box
0137NB	130 ml	Angio.	Mark III & Mark IV	Fill Tube	50-unit box
0137NN	150 ml	Angio.	Mark V, V+ & Provis	'Short nose' version + Fixed luer lock + Fill Tube	50-unit box
0137SC	150 ml	Angio.	Mark V, V+ & Provis	'Long nose' version + spinlock luer + Fill Tube	50-unit box
0137ND	200 ml	Angio.	Mark V, V+ & Provis	Fill Tube	50-unit box
0137NH	200 ml	Angio.	Mark V, V+ & Provis	-	50-unit box
<b>Covidien injectors (ex-Liebel Flarsheim)</b>					
0137NB	130 ml	CT	Optistat	Fill Tube	50-unit box
0137NE	150 ml	CT & Angio	Angiomat 6000	Fill Tube	50-unit box
0137NE	150 ml	CT & Angio.	Angiomat Illumena with pressure jacket 0700TB	Fill Tube	50-unit box
0137NA	100 ml	Angio.	Angiomat 3000 & Multijet	Fill Tube	50-unit box
0700TB	150 ml	CT & Angio	Pressure jacket for Angiomat Illumena		1-unit box
<b>Medtron injectors</b>					
0137NB	130 ml	CT, MRT & Angio	Injektron 82CT; CT2; MRT & 82HP	Fill Tube	50-unit box
0137ND	200 ml	CT, MRT & Angio	Injektron 82CT; CT2; MRT & 82HP	Fill Tube	50-unit box
0137NS	Kit 2 x 200 ml	CT	Injektron 82CT & CT2	Filling & injection set + 2 Fill tubes + 1 x 150 cm straight line	20-unit box
0137NT	Kit 130 ml & 200 ml	MRT	Injektron MRT	Filling & injection set + 2 Fill tubes + 1 x 150 cm straight line	20-unit box
<b>Nemoto injectors</b>					
0254NE	60 ml	CT	Dual Shot	Fill Tube	50-unit box
0254NF	100 ml	CT	Dual Shot, A300CT, A60CT	Fill Tube	50-unit box
0254NG	200 ml	CT	Dual Shot, A300CT, A60CT	Fill Tube	50-unit box
0254NA	Kit 60 ml & 200 ml	CT	Dual Shot, A300CT, A60CT	150 cm Y-connection line + 2 universal spikes	20-unit box
0254NB	Kit 60 ml	CT & MRT	CT > Dual Shot, A300CT, A60CT MRT > Sonic shot	150 cm Y-connection line + 1 universal spike	20-unit box
0254NC	Kit 200 ml	CT	Dual Shot, A300CT, A60CT	150 cm coiled connection line + 1 universal spike	20-unit box
0254ND	Kit 2 x 60 ml	CT & MRT	CT > Dual Shot, A300CT, A60CT MRT > Sonic shot	150 cm Y-connection line + 2 spikes (universal, mini)	20-unit box

Please contact **PEROUSE MEDICAL** for non sterile and OEM products.

**PEROUSE**  
MEDICAL

135, route Neuve - 69540 IRIGNY - FRANCE  
Phone +33 472 397 414 - Fax +33 478 518 967  
e-mail: [exportservice@sedat.com](mailto:exportservice@sedat.com) - [www.sedat.com](http://www.sedat.com)

